


**AF-1™ spheric**  
Hydrophobic IOL

Hydrophobic Acrylic IOL  
as low as 2.8 mm incision

**VA-60BB**



- ▶ **Product Configuration**  
Excellent centration and intraocular stability
- ▶ **Multiuse Injector System**  
Convenient lens delivery in HOYA's proprietary lens locking case
- ▶ **Wide Power Range**  
Available from -7.0 to +40.0 D



Multiuse Injector  
IS Blue Series Model ISH001




Cartridge C1

**HOYA**  
SURGICAL OPTICS

# AF-1<sup>TM</sup> spheric

Hydrophobic IOL

## VA-60BB



| Model                             | VA-60BB   |
|-----------------------------------|---|
| Specification                     | UV-Filtering <sup>1</sup>   |
| Optic Material                    | Hydrophobic Acrylic   |
| Optic Design                      | Spherical   |
| Manufacturing                     | Lathe-cut   |
| Haptic Material                   | PMMA chemically bonded  |
| Haptic Configuration              | Hopper loop, 5° angulation  |
| Dimensions (Optic/OAL)            | 6.0 mm / 12.5 mm  |
| Power                             | +4.0 to +40.0 D (0.5 D steps)<br>-7.0 to +3.0 D (1.0 D steps)   |
| A-Constant <sup>2</sup>           | 118.7   |
| IOL Master Constants <sup>3</sup> | Haigis a0 = 1.24 a1 = 0.4 a2 = 0.1<br>Hoffer Q pACD = 5.43<br>Holladay 1 sf = 1.65<br>SRK/T A = 118.7<br>SRK II A = 118.9 |
| Implantation System               | Multiuse Injector: IS Blue Series<br>Model Name: ISH001<br>Power: +4.0 to +26.0 D<br>Cartridge: C1                        |
| Incision                          | as low as 2.8 mm  |

- The AF-1 VA-60BB achieves a transmittance factor, which is close to that of the human lens. It blocks most ultraviolet light.
- This A-constant number is presented only as a guideline for lens power calculations. It is recommended that A-constant measurements be based on the surgeon's experience and measuring equipment.
- <http://ocusoft.de/ulib/c1.htm> (Stand: 26.04.2016)

HOYA Surgical Optics names and logos are registered trademarks of HOYA Surgical Optics, Inc. © 2016 HOYA Surgical Optics, Inc. All rights reserved.


HOYA Surgical Optics GmbH  
EMEA Headquarters  
De-Saint-Exupéry-Straße 8  
60549 Frankfurt am Main  
Germany

ifra-info@HOYA.com  
HOYA.com/SurgicalOptics


Singularly Focused. Globally Powered.<sup>TM</sup>

## Hydrophobic Acrylic IOL as low as 2.8 mm incision


Step A




Step B



Step C



Step D



Step A

Infuse OVD into the cartridge that sufficient amounts of OVD is inside of cartridge from tapered part to the end.

Step B

Fold the IOL by pushing the lever on the upper end of the provided lens case and pick up with forceps.

Step C

Hold the cartridge with the lens diagram facing upward, and then load the folded IOL into the cartridge.

Step D

Push the lens forward with leading haptic extending until the lens is placed in position as shown on the cartridge.



Insert the tip of the cartridge into the eye through the incision, keeping the bevel of the cartridge down.



Rotate the plunger clockwise and start releasing the lens into the eye. The leading haptic should be extending to the left (surgeon view).



After the leading haptic goes beneath the capsulorhexis, rotate the main body of the handpiece clockwise to release the IOL. The rotation direction is also indicated by an arrow on the main body.



Continuously rotate the plunger clockwise to completely release the optical portion of the lens from the cartridge and pull out the cartridge from the incision. The trailing haptic should be remained outside of the eye through the incision.



As carefully observing that the lens has recovered from its folded state in the eye without turning over, adjust the position of the lens to place the trailing haptic in the capsular bag, using a hook or other devices.

**HOYA**  
SURGICAL OPTICS