

REVERSO® MONOFOCAL

TECHNICAL SPECIFICATIONS



DESIGNATION	TECHNICAL SPECIFICATIONS
Lens type	For implantation into the ciliary sulcus, in pseudophakic patients
Optic diameter	6.50 mm
Overall diameter	13.80 mm
Design	One piece round edge on 360°
Optic design	Spherical, convex anterior surface Concave posterior surface
Angulation	10°
Material	25% hydrophilic acrylic
Dioptric powers	On request: from -6.0D to +6.0D by 0.5D
Refractive index	1,46
Sterilization	Steam
Recommended incision size	< 2.0 mm

MADE
IN FRANCE
CE marked

► Reverso calculator available on our website www.cristalens.fr

Hydrophilic Acrylic

Correction of a potential refractive error in pseudophakic patients

Micro incision

Monofocal

Reference: REVERSO