

MSP™

A COMPLETE LINE OF

OPHTHALMIC KNIVES



MSP™ Ophthalmic Microsurgical Knives

- Trusted In Ophthalmology
- Proprietary Finishing Technology
- Precision Manufacturing

SURGICAL
SPECIALTIES
CORPORATION™

Round Stock

Stab Incision

Stab knives and blades for initial stab incisions. 6 units per box.



Knife	Description
7503	15° Blade
7509	22.5° Blade
7505	30° Blade

Blade	Description
7502	15° Blade

Microsurgical Controlled Depth

Precise microsurgical knives and blades for controlled depth incisions. 6 units per box.



7514 15°, 3.0mm

7513 15°, 3.0mm



7516 15°, 5.0mm

7515 15°, 5.0mm



7523 22.5°, 4.0mm

7522 22.5°, 4.0mm



7531 30°, 3.5mm

Crescent – Scleral Tunnel

6 units per box.



66WUDK Angled, Bevel Up

Phaco – Slit

Precise width phaco-slit knives are designed to facilitate optimal fit for phaco tips. 6 units per box.



Knife

Description

56B220BDK	2.2mm, Angled, Double Bevel
56B240BDK	2.4mm, Angled, Double Bevel
56B250RDK	2.5mm, Angled, Bevel Up
56B275RDK	2.75mm, Angled, Bevel Up
56B280RDK	2.8mm, Angled, Bevel Up
56B300RDK	3.0mm, Angled, Bevel Up
56B320RDK	3.2mm, Angled, Bevel Up
56B350RDK	3.5mm, Angled, Bevel Up

Implant

Implant knives allow for a standard, reproducible widening of the incision to accommodate IOL width. 6 units per box.



56B350RTDK	3.5mm, Round Tip, Angled, Bevel Up
56B560RTDK	5.6mm, Round Tip, Angled, Bevel Up

MVR

6 units per box



5546	19 Gauge, Straight
5547	20 Gauge, Straight
5547B	20 Gauge, Angled

Flat Stock

Scleral

6 units per box.



Knife	Description
57DK	57 Scleral



Blade	Description
57	57 Scleral

Spoon

6 units per box.



SF280BDK	Medium, Angled
----------	----------------



SF150	Small, Straight
-------	-----------------

Mini-Edge

144 units. (12 boxes of 12 units per box)



64	Round Tip
----	-----------



65	Lance Tip
----	-----------



67	Curved Tip
----	------------



69	Sharp All Around
----	------------------

Needle

6 units per box.



59S	3.0mm
-----	-------

MSP™

Flat Stock

Dissection

6 units per box.



Blade	Description
52S	Small
52SR	Small, Special Grind

Lamellar

6 units per box.



6660	Angled, 6600 Blade
------	--------------------

Stainless Steel Handles

Sold individually.



3KM	7.5cm, Round Knurled
-----	----------------------



3K	10cm, Round Knurled
----	---------------------



3H	10cm, Hex Handle
----	------------------



5S	3.5 inch, Self Locking
----	------------------------



Visit SurgicalSpecialties.com For More Information

Surgical Specialties Corporation, 50 Braintree Hill Office Park, Suite 101, Braintree, MA 02184, U.S.A.
Phone: 610.404.3520 • 877.991.1110 (USA and Canada only) • Fax: 610.404.4010
All trademarks are owned by the company and its affiliates. ©2001-2015 Surgical Specialties Corporation.
All Rights Reserved. SSC-046 R1 10/16

**SURGICAL
SPECIALTIES**
CORPORATION™