

INTRAOCULAR LENS


**ARTIS<sup>®</sup> PLE**

**HYDROPHOBIC ACRYLIC**



Micro-incision

*Preloaded*



Made in France

CE 0459



DESIGNATION	TECHNICAL SPECIFICATIONS
<b>Lens Type</b>	For implantation in the capsular bag
<b>Optic diameter</b>	6.15 mm (from oD to +9.5D) 6.00 mm (from +10D to +25D) 5.80 mm (from +25.5D to +35D)
<b>Overall diameter</b>	11.00 mm (from oD to +9.5D) 10.79 mm (from +10D to +25D) 10.50 mm (from +25.5D to +35D)
<b>Design</b>	One piece square edge on 360°
<b>Optic design</b>	Aspherical on the posterior surface Compensation of corneal aberrations
<b>Angulation</b>	5°
<b>Material</b>	Hydrophobic acrylic for micro-incision
<b>Dioptric powers</b>	From oD to +35D by 0.5D
<b>Estimated A-Constant (SRK-T)</b>	119.3 Ultrasound biometry 119.7 Interference laser biometry
<b>Suggested Anterior Chamber Depth (ACD)</b>	5.77 mm Ultrasound biometry 6.03 mm Interference laser biometry
<b>Refractive index</b>	1.54
<b>Sterilization</b>	Ethylene oxide
<b>Recommended incision size</b>	< 2.0 mm oval shape

04/05/2015 DOC\_10\_PHOBIC



Ophthalmologic  
surgery  
**Cristalens**

12 allée Rosa Luxemburg - BP 50240 Eragny  
95615 Cergy Pontoise cedex - France

Tél : + 33 (0) 1 34 32 30 50 - Fax : + 33 (0) 1 34 32 30 60



**Cristalens Industrie**

4 rue Louis de Broglie - 22300 Lannion - France

Tél : + 33 (0) 2 96 48 92 92 - Fax : + 33 (0) 2 96 48 97 87

[www.cristalens.fr](http://www.cristalens.fr)

Cristalens